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Responsabile: Prof. Andrea Marzoli

Tecnici: Dr. Daria Pasqual

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X-RAY FLUORESCENCE (XRF) SPECTROMETRY LABORATORY

The instrument is a WDS sequential Philips PW2400 spectrometer, operating under vacuum conditions, using the SuperQ software, and equipped with:

- 3 kW Rh X-ray tube
- 5 analyzing crystals (LiF220, LiF200, Ge, PE, TIAP)
- 2 detectors (flow counter and scintillator)
- 3 collimators (150 μm , 300 μm and 700 μm)
- 4 filters (Al 200 μm , Brass 100 μm , Pb 1000 μm and Brass 300 μm)
- sample changer Philips PW2510 with 30 sample holders.

Routine analysis of a sample

- e. Determination of loss on ignition (L.O.I).
- f. Use of the calcined sample to prepare the glass bead.
- g. Sample analysis using the WDS sequential Philips PW2400 spectrometer.
- h. Determination of Fe^{2+} content via permanganometry.

a. Determination of loss on ignition (L.O.I): a known amount of sample is heated in a furnace at 860°C for about 20 mins, then at 980°C for about 2 hours.

b. Sample preparation for XRF analysis: samples for XRF analysis are glass beads prepared using calcined sample diluted with flux di-lithium tetraborate $\text{Li}_2\text{B}_4\text{O}_7$ (1:10 ratio) and melted with a fluxer Claisse Fluxy (~1150°C).

c. Sample analysis: Chemical elements analyzed are Si, Ti, Al, Fe, Mn, Mg, Ca, Na, K, P, Sc, V, Cr, Co, Ni, Cu, Zn, Ga, Rb, Sr, Y, Zr, Nb, Ba, La, Ce, Nd, Pb, Th and U.

Reference standards are natural geological samples (Geostandards Newsletter, Vol. XVIII, Special Issue, July 1994, K. Govindaraju, ed.).

Instrumental precision (defined by several measurements performed on the same sample) is within 0.6% relative for major elements, and within 3% relative for trace elements.

Detection limits are within 0.01% for major elements Al, Mg and Na, within 0.2% for Si, within 0.005% for Ti, Fe, Mn, Ca, K and P. Limits for trace elements are (in ppm): Sc=5, V=5, Cr=6, Co=3, Ni=3, Cu=3, Zn=3, Ga=5, Rb=3, Sr=3, Y=3, Zr=3, Nb=3, Ba=10, La=10, Ce=10, Nd=10, Pb=5, Th=3, U=3.

d. Determination of Fe^{2+} content: Fe measured with XRF is expressed as Fe_2O_3 , while the FeO content is determined through potassium permanganate titration (permanganometry, Pratt method).

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